

Product Datasheet

NF2/Merlin antibody (orb640943)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to NF2/Merlin |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IP, WB |
| Immunogen | Recombinant full length Human NF2/Merlin. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 80kDa |
| Uniprot ID | P35240 |
| Dilution Range | WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200 |
| Expiration Date | 12 months from date of receipt. |