

Product Datasheet

Sam 68 antibody (orb641804)

Description	Rabbit polyclonal antibody to Sam 68
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human KHDRBS1/Sam68(NP_006550.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	68kDa
Uniprot ID	Q07666
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:100 - 1:200, IF1:50 - 1:200
Expiration Date	12 months from date of receipt.