

Product Datasheet

Contactin-2 antibody (orb642022)

Description	Rabbit polyclonal antibody to Contactin-2
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human CNTN2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 130 kDa
Uniprot ID	Q02246
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.