

Product Datasheet

PDGFR- alpha (phospho-Y754) antibody (orb642237)

Description	Rabbit polyclonal antibody to PDGFR-alpha (phospho-Y754)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y754 of human PDGFRA.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 180 kDa
Uniprot ID	P16234
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.