

Product Datasheet

HA-tag antibody (orb642261)

Description	Rabbit polyclonal antibody to HA-tag
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide HA-tag.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:2000-1:10000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.