

Product Datasheet

HM74 (T327) antibody (orb642292)

Description	Rabbit polyclonal antibody to HM74 (T327)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 300-350 of Human HM74.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 45 kDa
Uniprot ID	Q8TDS4 , P49019
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.