

Product Datasheet

TK1 (phospho-Ser13) antibody (orb642302)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TK1 (phospho-Ser13) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic phosphopeptide derived from human TK1 around the phosphorylation site of Serine 13. |
| Preservatives | 1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 26 kDa |
| Uniprot ID | P04183 |
| Dilution Range | WB: 1:500-1:1000 |
| Expiration Date | 12 months from date of receipt. |