

Product Datasheet

Cdk9 (phospho-T186) antibody (orb642317)

Description	Rabbit polyclonal antibody to Cdk9 (phospho-T186)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Cdk9 around the phosphorylation site of Threonine 186
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 42, 55 kDa
Uniprot ID	P50750
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.