

## Product Datasheet

# NMDA Epsilon 2 (phospho-Y1336) antibody (orb642320)

<b>Description</b>	Rabbit polyclonal antibody to NMDA Epsilon 2 (phospho-Y1336)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Synthetic phosphopeptide derived from human NMDA $\epsilon$ 2 around the phosphorylation site of Tyrosine1336.
<b>Preservatives</b>	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 190 kDa
<b>Uniprot ID</b>	<a href="#">Q13224</a>
<b>Dilution Range</b>	WB: 1:500-1:1000
<b>Expiration Date</b>	12 months from date of receipt.