

Product Datasheet

SREBP-1 (phospho-S439) antibody (orb642321)

Description	Rabbit polyclonal antibody to SREBP-1 (phospho-S439)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human SREBP-1 around the phosphorylation site of Serine 439
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 121 kDa
Uniprot ID	P36956
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.