

Product Datasheet

Myc-Tag antibody (orb642331)

Description	Rabbit polyclonal antibody to Myc-Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide Myc-Tag.
Preservatives	1 mg/ml in Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	N/A
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.