

Product Datasheet

PKC Zeta (phospho-T560) antibody (orb642354)

Description	Rabbit polyclonal antibody to PKC Zeta (phospho-T560)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human PKC ζ around the phosphorylation site of Threonine 560.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 70 to 85 kDa
Uniprot ID	Q05513
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:10-1:100
Expiration Date	12 months from date of receipt.