

Product Datasheet

His-tag (4C2) antibody (orb642375)

Description	Mouse monoclonal antibody to His-tag (4C2)
Species/Host	Mouse
Reactivity	All
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Synthetic peptide His-tag.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB:1:5000-10000IP 1:50 - 1:100
Expiration Date	12 months from date of receipt.