

Product Datasheet

HSP27 (phospho-S15) antibody (orb642394)

Description	Rabbit polyclonal antibody to HSP27 (phospho-S15)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human HSP27 around the phosphorylation site of Serine 15.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 27 kDa
Uniprot ID	P04792
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:10-1:100
Expiration Date	12 months from date of receipt.