

Product Datasheet

PKC alpha (phospho-T638) antibody (orb642400)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to PKC alpha (phospho-T638) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, IP, WB |
| Immunogen | Synthetic phosphopeptide derived from human PKC α around the phosphorylation site of Threonine 638. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 80 kDa |
| Uniprot ID | P17252 |
| Dilution Range | WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:10-1:100 |
| Expiration Date | 12 months from date of receipt. |