

Product Datasheet

Smad3 (phospho-S213) antibody (orb642407)

Description	Rabbit polyclonal antibody to Smad3 (phospho-S213)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic phosphopeptide derived from human Smad3 around the phosphorylation site of Serine 213.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 48, 55 kDa
Uniprot ID	P84022
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.