

Product Datasheet

TIMP-1 (G94) antibody (orb643063)

Description	Rabbit polyclonal antibody to TIMP-1 (G94)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to amino acids 61-110 of Human TIMP-1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 23 kDa
Uniprot ID	P01033
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.