

Product Datasheet

TCF-4/12 (S616) antibody (orb643218)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to TCF-4/12 (S616) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Synthetic peptide, corresponding to amino acids 580-630 of Human TCF-4. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 60, 72 kDa |
| Uniprot ID | P15884, Q99081 |
| Dilution Range | WB: 1:500-1:1000, IHC-P: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |