

Product Datasheet

AMPK Gamma 1/2/3 (K78) antibody (orb643586)

Description	Rabbit polyclonal antibody to AMPK Gamma 1/2/3 (K78)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic peptide, corresponding to amino acids 50-100 of Human AMPK γ 1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 38, 54, 63 kDa
Uniprot ID	P54619 , Q9UGI9 , Q9UGJ0
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.