

Product Datasheet

Caspase 2 (K152) antibody (orb643843)

Description	Rabbit polyclonal antibody to Caspase 2 (K152)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 140-180 of Human Caspase 2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 51 kDa
Uniprot ID	P42575
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.