

Product Datasheet

eEF2K (phospho-S366) antibody (orb644397)

Description	Rabbit polyclonal antibody to eEF2K (phospho-S366)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human eEF2K around the phosphorylation site of Serine 366.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 105 kDa
Uniprot ID	O00418
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:10-1:100
Expiration Date	12 months from date of receipt.