

Product Datasheet

ICAM-1 (phospho-Y512) antibody (orb644424)

Description	Rabbit polyclonal antibody to ICAM-1 (phospho-Y512)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human ICAM-1 around the phosphorylation site of Tyrosine 512.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 58, 96 kDa
Uniprot ID	P05362
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.