

## Product Datasheet

# Integrin beta 3 (phospho-Y773) antibody (orb644429)

<b>Description</b>	Rabbit polyclonal antibody to Integrin beta 3 (phospho-Y773)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC
<b>Immunogen</b>	Synthetic phosphopeptide derived from human Integrin $\beta$ 3 around the phosphorylation site of Tyrosine 773.
<b>Preservatives</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 87, 110 kDa
<b>Uniprot ID</b>	<b>P05106</b>
<b>Dilution Range</b>	IHC-P: 1:50-1:200
<b>Expiration Date</b>	12 months from date of receipt.