

## Product Datasheet

### MEF-2 (phospho-T312) antibody (orb644441)

<b>Description</b>	Rabbit polyclonal antibody to MEF-2 (phospho-T312)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, IP
<b>Immunogen</b>	Synthetic phosphopeptide derived from human MEF-2 around the phosphorylation site of Threonine 312.
<b>Preservatives</b>	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 55 kDa
<b>Uniprot ID</b>	<a href="#">Q02078</a>
<b>Dilution Range</b>	IHC-P: 1:50-1:200IP: 1:10-1:100
<b>Expiration Date</b>	12 months from date of receipt.