

Product Datasheet

VEGFR2 (phospho-Y1214) antibody (orb644501)

Description	Rabbit polyclonal antibody to VEGFR2 (phospho-Y1214)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic phosphopeptide derived from human VEGFR2 around the phosphorylation site of Tyrosine 1214.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 230 kDa
Uniprot ID	P35968
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.