

Product Datasheet

E2F1 (phospho-T433) antibody (orb644519)

Description	Rabbit polyclonal antibody to E2F1 (phospho-T433)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human E2F-1 around the phosphorylation site of Threonine 433.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 47 kDa
Uniprot ID	Q01094
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.