

Product Datasheet

p70 S6K (phospho-S418) antibody (orb644549)

Description	Rabbit polyclonal antibody to p70 S6K (phospho-S418)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human p70 S6K around the phosphorylation site of Serine 418.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 59,70,85 kDa
Uniprot ID	P23443
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.