

Product Datasheet

Troponin I3 (phospho-S22/S23) antibody (orb644553)

Description	Rabbit polyclonal antibody to Troponin I3 (phospho-S22/S23)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human Troponin I3 around the phosphorylation site of Ser22/23.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 28 kDa
Uniprot ID	P19429
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.