

Product Datasheet

Ets-1 (phospho-T38) antibody (orb644567)

Description	Rabbit polyclonal antibody to Ets-1 (phospho-T38)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic phosphopeptide derived from human ETS1 around the phosphorylation site of Threonine 38.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 54 kDa
Uniprot ID	P14921
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.