

Product Datasheet

LAT (phospho-Y191) antibody (orb644573)

Description	Rabbit polyclonal antibody to LAT (phospho-Y191)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human LAT (Isoform Short) around the phosphorylation site of Tyrosine 191.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 38 kDa
Uniprot ID	O43561
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.