

Product Datasheet

Cdk5 (phospho-Y15) antibody (orb644624)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Cdk5 (phospho-Y15) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic phosphopeptide derived from human Cdk5 around the phosphorylation site of Tyrosine15. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 33 kDa |
| Uniprot ID | Q00535 |
| Dilution Range | WB: 1:500-1:1000 |
| Expiration Date | 12 months from date of receipt. |