

Product Datasheet

Dok-2 (phospho-Y299) antibody (orb644716)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to Dok-2 (phospho-Y299) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | Synthetic phosphopeptide derived from human Dok-2 around the phosphorylation site of Tyrosine 299. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 56 kDa |
| Uniprot ID | O60496 |
| Dilution Range | WB: 1:500-1:1000, IHC-P: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |