

Product Datasheet

Shb (phospho-Y246) antibody (orb644732)

Description	Rabbit polyclonal antibody to Shb (phospho-Y246)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human Shb around the phosphorylation site of Tyrosine 246.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~55kDa
Uniprot ID	Q15464
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.