

Product Datasheet

GATA-1 (phospho-S142) antibody (orb644739)

Description	Rabbit polyclonal antibody to GATA-1 (phospho-S142)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human GATA1 around the phosphorylation site of Serine 142.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 43 kDa
Uniprot ID	P15976
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200IP: 1:50-1:200
Expiration Date	12 months from date of receipt.