

Product Datasheet

Smad1 (phospho-S187) antibody (orb644764)

Description	Rabbit polyclonal antibody to Smad1 (phospho-S187)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic phosphopeptide derived from human Smad1 around the phosphorylation site of Serine 187.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	Q15797
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:200 ICC: 1:50-1:200 IP: 1:50-1:200
Expiration Date	12 months from date of receipt.