

Product Datasheet

Catenin- beta (phospho-S33/S37/T41) antibody (orb644797)

Description	Rabbit polyclonal antibody to Catenin- beta (phospho-S33/S37/T41)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A phospho specific peptide corresponding to residues surrounding Serine 33 and 37 and threonine 41 of human Catenin- β .
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 92 kDa
Uniprot ID	P35222
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.