

Product Datasheet

PKC-pan antibody (orb646345)

Description	Rabbit polyclonal antibody to PKC-pan
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to Human PRKCA.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 82 kDa
Uniprot ID	Q02156 , P17252 , P05771 , Q05513 , Q05655 , Q04759 , P05129 , P24723
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.