

Product Datasheet

PSMA3 antibody (orb646423)

Description	Rabbit polyclonal antibody to PSMA3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant full length Human PSMA3.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 28 kDa
Uniprot ID	P25788
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200IP 1:50 - 1:100
Expiration Date	12 months from date of receipt.