

Product Datasheet

ASK1 (phospho-S1033) antibody (orb646452)

Description	Rabbit polyclonal antibody to ASK1 (phospho-S1033)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S1033 of human MAP3K5(NP_005914.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	155kDa
Uniprot ID	Q99683
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.