

## Product Datasheet

### EIF4G1 (phospho-S1148) antibody (orb646477)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to EIF4G1 (phospho-S1148)  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | A synthetic peptide corresponding to residues in Human EIF4G1 around the phosphorylation site of S1148.                                 |
| <b>Preservatives</b>       | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2   |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | ~ 175 kDa   |
| <b>Uniprot ID</b>          | <a href="#">Q04637</a>  |
| <b>Dilution Range</b>      | WB: 1:500-1:1000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |