

Product Datasheet

DAPK1 (phospho-S736) antibody (orb646503)

Description	Rabbit polyclonal antibody to DAPK1 (phospho-S736)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic phosphopeptide derived from human DAPK1 around the phosphorylation site of Ser736.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 160 kDa
Uniprot ID	P53355
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.