

Product Datasheet

Cleaved-Notch 1 (V1754) antibody (orb646784)

Description	Rabbit polyclonal antibody to Cleaved-Notch 1 (V1754)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic peptide, corresponding to amino acids 1700-1800 of Human Notch 1.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 110 kDa
Uniprot ID	P46531
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.