

Product Datasheet

DNA pol Delta cat antibody (orb647134)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to DNA pol Delta cat |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Immunogen | Recombinant full length Human DNA pol δ cat. |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | ~ 125 kDa |
| Uniprot ID | P28340 |
| Dilution Range | WB: 1:500-1:2000, IHC-P: -P/, IF: 1:50-1:200 |
| Expiration Date | 12 months from date of receipt. |