



Product Datasheet

CGRP antibody (orb647197)


Description	Rabbit polyclonal antibody to CGRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC
Immunogen	Recombinant protein of human CALCB.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 14 kDa
Uniprot ID	P10092
Dilution Range	WB 1:500 - 1:2000 , IHC-P: -P 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin embedde...



Immunohistochemistry of paraffin embedde...



Western blot analysis of extract of vari...