

Product Datasheet

GABAA R beta 2 antibody (orb647280)

Description	Rabbit polyclonal antibody to GABAA R beta 2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length Human GABAA R β 2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 59 kDa
Uniprot ID	P47870
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.