

Product Datasheet

EGLN1 antibody (orb648207)

Description	Rabbit polyclonal antibody to EGLN1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human EGLN1/EGLN2(NP_071334.1/NP_444274.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	44kDa
Uniprot ID	Q96KS0 , Q9GZT9
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200
Expiration Date	12 months from date of receipt.