

Product Datasheet

Fyn (phospho-S21) antibody (orb648588)

Description	Rabbit polyclonal antibody to Fyn (phospho-S21)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S21 of human Phospho-Fyn-S21(NP_002028.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	60kDa
Uniprot ID	P06241
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.