

Product Datasheet

NDRG1 (phospho-S330) antibody (orb648598)

Description	Rabbit polyclonal antibody to NDRG1 (phospho-S330)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S330 of human NDRG1(NP_001128714.1).
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q92597
Dilution Range	WB1:500 - 1:2000
Expiration Date	12 months from date of receipt.