

Product Datasheet

NEDD4L (phospho-S448) antibody (orb648601)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to NEDD4L (phospho-S448) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | A synthetic phosphorylated peptide around S448 of human NEDD4L |
| Preservatives | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Form/Appearance | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2 |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | 130kDa |
| Uniprot ID | Q96PU5 |
| Dilution Range | WB1:500 - 1:2000 |
| Expiration Date | 12 months from date of receipt. |