

Product Datasheet

Histone H3 (aDi-Methyl R2) antibody (orb648645)

Description	Rabbit polyclonal antibody to Histone H3 (aDi-Methyl R2)
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IP, WB
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1-100 of Human Histone H3, aDi-methylated at R2.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~15 kDa
Uniprot ID	Q16695
Dilution Range	WB1:500 - 1:2000, IHC-P: -P1:50 - 1:200, IF1:50 - 1:200IP1:50 - 1:200, ChIP1:20 - 1:100, ChIPseq1:20 - 1:100
Expiration Date	12 months from date of receipt.